

# The First Original is the Best!

The Two-in-one (RO&UP combined) system, which produces RO and UP water in one system, was firstly developed in 1993 by the **Human Corporation**.



Not only was this a world first but also it set the standard for other manufacturers to follow with its compact design and effective cost. Here is one of the systems called "P.NIX power" which is fully CE certified, manufactured under strictly controlled ISO 9001 standards.

Standard features include self-test, various kinds of auto cleaning and recycling functions by timer as well as by the water quality. To maintain optimum water quality is ensured with automatic sensor calibration and 2 steps of filter exchange indicators.

It is further enhanced with programmable timed dispensing modes and computer interfaced monitoring systems on critical stages of the purification process.